

PROTOCOL FOR STS-PCRSEQ

Strand-and-Transcript-Specific (STS)-PCRSeq is a method to assess the efficiency of RNA editing in a discrete set of transcripts. Strand-specific primers are used to generate cDNA. The cDNA is then amplified by PCR and subjected to high-throughput sequencing. While the method has been used for assessing Arabidopsis organelle editing extent, the method could be applied to any biological system in which editing of multiple transcripts occurs.

Several steps are necessary for using STS-PCRseq to measure plant organelle editing extent

1-Extract RNA from plants.

Best method to produce high quality RNA consists in extracting RNA with Trizol then follow by using the Pure Link RNA Mini Kit (Invitrogen).

1. Grind 2-3 young leaves in liquid nitrogen, re-suspend in 800 μ l of Trizol.
2. Add 160 μ l of chloroform, centrifuge 13000 rpm in a microfuge for 15 minutes (cold room)
3. Recover upper phase (~400 μ l)
4. Then follow the following protocol from the Pure Link RNA Mini Kit (Invitrogen).
 - 4.1. Add one volume 70% ethanol (~400 μ l) to the recovered upper phase
 - 4.2. Vortex to mix thoroughly and to disperse any visible precipitate that may form after adding ethanol.
 - 4.3. Transfer up to 700 μ l of the sample (including any remaining precipitate) to the Spin Cartridge (with the Collection Tube).
 - 4.4. Centrifuge at 12,000 \times g for 14 seconds at room temperature. Discard the flow-through, and reinsert the Spin Cartridge into the same Collection Tube.
 - 4.5. Repeat Steps 4.3–4.4 until the entire sample has been processed.
 - 4.6. Add 700 μ l Wash Buffer I to the Spin Cartridge. Centrifuge at 12,000 \times g for 15 seconds at room temperature. Discard the flow-through and the Collection Tube. Place the Spin Cartridge into a new Collection Tube.
 - 4.7. Add 400 μ l Wash Buffer II with ethanol (page 11) to the Spin Cartridge.
 - 4.8. Centrifuge at 12,000 \times g for 15 seconds at room temperature. Discard the flow-through and reinsert the Spin Cartridge into the same Collection Tube.
 - 4.9. Repeat Steps 4.7–4.8 once.
 - 4.10. Centrifuge the Spin Cartridge at 12,000 \times g for 1-2 minutes to dry the membrane with attached the RNA. Discard the Collection Tube and insert the Spin Cartridge into a Recovery Tube.

- 4.11. Add 30 μ L RNase-Free Water to the center of the Spin Cartridge (see Elution Parameters, page 12).
- 4.12. Incubate at room temperature for 1 minute.
- 4.13. Centrifuge the Spin Cartridge for 2 minutes at $\geq 12,000 \times g$ at room temperature to elute the RNA from the membrane into the Recovery tube.
5. DNase the RNA with Turbo DNA free kit from Ambion
 - 5.1. Add 3 μ L (0.1 volume) 10 \times TURBO DNase Buffer and 1 μ L TURBO DNase to the RNA, and mix gently
 - 5.2. Incubate at 37°C for 30 min
 - 5.3. Add 1 μ L TURBO DNase to the RNA, and mix gently
 - 5.4. Incubate at 37°C for 30 min
 - 5.5. Add 3 μ L resuspended DNase Inactivation Reagent and mix well
 - 5.6. Incubate 5 min at room temperature, mixing occasionally.
 - 5.7. Centrifuge at 10,000 $\times g$ for 1.5 min and transfer the RNA to a fresh tube

2- Reverse-transcribe organelle transcripts with strand-specific primers

Best results with the Superscript III Reverse Transcriptase (Invitrogen)

1. Mix all the plastid reverse primers in one tube (2 μ M) and all the mitochondrial reverse primers in another tube (2 μ M) (all the primers are given in the Table at the end of the protocol)
2. Perform the RT reaction with the Superscript III Reverse Transcriptase (Invitrogen)
 - 2.1. Add the following components to a nuclease-free microcentrifuge tube:
 - 1 μ L of either plastid or mitochondrial gene-specific primer mix
 - 1 μ L total DNAsed RNA
 - 1 μ L 10 mM dNTP Mix (10 mM each dATP, dGTP, dCTP and dTTP at neutral pH)
 - Sterile, distilled water to 13.5 μ L
 - 2.2. Heat mixture to 65°C for 5 minutes and incubate on ice for at least 1 minute
 - 2.3. Collect the contents of the tube by brief centrifugation and add:
 - 4 μ L 5X First-Strand Buffer
 - 1 μ L 0.1 M DTT
 - 0.5 μ L RNaseOUT™ Recombinant RNase Inhibitor (Cat. no. 10777-019,40 units/ μ L).

1 μ l of SuperScript™ III RT (200 units/ μ l)*

*** I generally use RT diluted 20 times, which is fine to reverse transcribe organelle transcripts and saves a lot of money.**


2.4. Mix by pipetting gently up and down.

2.5. Incubate at 50°C for 60 minutes.

2.6. Inactivate the reaction by heating at 70°C for 15 minutes.

3- Amplification of all mitochondrial and plastid cDNAs

Originally, all the PCR reactions were performed individually, one PCR reaction for one gene. Because there are 61 organelle genes (36 mitochondrial + 25 plastid) it means 61 PCR reactions. I recently developed a multiplex PCR strategy where I combine from 2 to 4 primer pairs, reducing the total number of PCR reactions to 25 (Primer combinations are given in a Table at the end of protocol). Significant amount of troubleshooting was necessary in order to find the appropriate primer combinations allowing amplifications to occur at roughly the same molar yield. An example is given below:



Negative of an Ethidium Bromide gel (1.5% Synergel + 0.7% Agarose)

Because of the constraints on the molar yield of the amplicons, the PCR primer pairs were chosen so that the resulting PCR products sizes were in a close range (see above gel). I use the Synergel agarose additive (Diversified Biotech) that allows improved separation and definition of DNA fragments.

1. Amplification of RT-PCR products. I use the Biomix™ Red (Bioline) which is a complete ready-to-use 2x reaction mix containing an ultra-stable Taq DNA polymerase. It contains an additional inert red dye that permits easy visualization and direct loading onto a gel.

1.1. PCR mix

BioMix Red	12.5 μ l
Primer combination	0.1 μ M
cDNA	1 μ l
Water (nuclease-free)	up to 25 μ l

1.2. Thermal cycling conditions

Step	Temperature °C	Time	Number of cycles
Initial denaturation	94	2 min	1
Denaturation	94	30s	45
Annealing	50	30s	
Extension	72	1 min/kb	
Final extension	72	10 min	1

2. Preparing a 1.5% Synergel + 0.7% Agrose gel


For a big tray gel, I need 300 ml of 1.5% Synergel + 0.7% Agrose or 4.5 g of Synergel and 2.1g of regular agarose

- 2.1. Add the Synergel™ powder to the beaker or flask and thoroughly intermix with the agarose.
- 2.2. Thoroughly saturate the mixed powders with ethanol. (Use enough ethanol so that the resulting slurry can just move about freely at the bottom of the container).
- 2.3. Add the required amount of running buffer: 300 ml of 1xTAE
- 2.4. Dissolve the suspension in a microwave
- 2.5. Pour the gel and set the combs.

4- Quantification of the RT-PCR products

The STS-PCRseq is aimed at measuring editing extent for every organelle gene transcript and untranslated regions that have been reported to contain edited sites. In order to achieve a similar accuracy for every edited site we mix the RT-PCR products in equimolar ratio. The quantification of RT-PCR products is done by using ImageJ software on the captured image of the Ethidium Bromide loaded gel such as the one shown above (generally I load 5 µl of the PCR reaction saving 20 µl for further use in the final mix). ImageJ is a public domain Java image processing program inspired by NIH Image for the Macintosh that can calculate area and pixel value statistics of user-defined selections (<http://imagej.nih.gov/ij/docs/intro.html>).

As an example I am showing the resulting image from the plot lane function obtained with one of the lane from the gel above and the three molecular ladder lanes



The number of pixels under the molecular ladder peaks allows the establishment of a linear relationship between pixel intensity and amount in ng. This relationship is then used to determine the amount of the RT-PCR products.

After quantifying all the RT-PCR products, they are mixed up in equimolar ratio (100 fmoles/RT-PCR)

5- Cleaning of the RT-PCR products mix

In my most recent STS-PCRseq experiment, the volume of the RT-PCR mixes ranged from 50 μ l to 200 μ l. Before sharing the mix by sonication, they have to be cleaned from the primers, nucleotides, salts contained in the PCR reactions. I use the PureLink™ quick PCR purification kit (Invitrogen).

- 5.1. Add 4 volumes of Binding Buffer B2 with isopropanol to 1 volume of a PCR sample. Mix well.
- 5.2. Load. Pipet the sample into a PureLink® Spin Column in a Collection Tube. Centrifuge the column at $>10,000 \times g$ for 1 minute. Discard the flow-through.
- 5.3. Wash. Re-insert the column into the Collection Tube and add 650 μ L Wash Buffer (W1) with ethanol. Centrifuge the column at $>10,000 \times g$ for 1 minute. Discard the flow-through and place the column in the same Collection Tube. Centrifuge the column at maximum speed for 2–3 minutes.
- 5.4. Elute. Place the column into a clean 1.7-mL Elution Tube (supplied with the kit). Add 50 μ L nuclease free-water to the center of the column. Incubate the column at room temperature for 1 minute. Centrifuge the column at maximum speed for 2 minutes. The elution tube contains the purified PCR product. Store the purified DNA at 4°C for immediate use or at -20°C for long-term storage.

6- Quantification of the RT-PCR products mix by Qubit dsDNA HS Assay Kit (molecular probes)

The Qubit assay is very sensitive and target-specific using fluorescent dyes that bind specifically to the target molecule, thus minimizing the effect of contaminants-including DNA or RNA-on the result. It requires the assay to be read in the Qubit fluorometer.





Figure from ThermoFisher Scientific website

After Qubit quantification, 1200ng of each RT-PCR mix was adjusted to a sample with a volume of 60 µl and a concentration of 20ng/ µl

7- Shearing of 1.1 µg of RT-PCR mix (55 µl) by sonication with a Covaris M220 Focused Ultrasonicator

5 µl of each RT-PCR mix is set aside to monitor the efficiency of the shearing. 55 µl of the remaining RT-PCR mix is put in a Covaris microTUBE™ and undergoes shearing according to the manufacturer recommendations to generate 300 bp fragments. After shearing of the mix I load a gel with neighboring lanes containing 5 µl of the mix un-sheared and 5 µl of the mix sheared



B: before shearing

A: after shearing

8- Preparation of a TruSeq DNA PCR-Free library

1 µg of sheared RT-PCR mix enters the workflow of the preparation of a TruSeqDNA PCR-Free library from the HT kit (Illumina).

Citation of this protocol. The reads processing and alignment have been described in the following paper: “Comprehensive High-Resolution Analysis of the Role of an Arabidopsis Gene Family in RNA Editing”. PLOS Genetics 9: e1003584. S Bentolila S, J Oh , MR Hanson, R Bukowsky R (2013). <http://www.ncbi.nlm.nih.gov/pubmed/23818871>.

Primers used for organelle RT-PCR.							
purpose	organelle	gene		primers			product size (bp)
editing	plastid	<i>accD</i>	F2	ATGGAAAAATCGTGGTTCAAATTTATGTT	R2	GTTTGTCTAGTCTAATTGAACCTCCCC	1638
editing	plastid	<i>atpF</i>	F2	TTAACCGATTCTTTCGTTCACITG	R2	CAAAACATCCCAATAATTTGCAATTA	520
editing	plastid	<i>clp</i>	F2	TATTGGCGTTCAAAAAGTACCTT	R2	GAAACCGCTACAAGAACAACAATTC	581
editing	plastid	<i>matK</i>	F3	ATAAAATTTCAAGGATATTTAGAGTTCG	R3	TTGACCAAATCAATTAAGATAAAGAATA	1501
editing	plastid	<i>ndhB</i>	F3	GCA TGTA CAGAA TGAAAA TTCA TTCTC	R2	AA TCGCAA TAA TCGGGTTCATT	1504
editing	plastid	<i>ndhD</i>	F	AA CA A CT C G A A G T A T G G G T C	R2	CTAA TGAGAGCACA AAA ATCA G G A T	1637
editing	plastid	<i>ndhF</i>	F	GGATCA TACCTTTCATTCCACTTC	R	GCA GCATGTATAAGAGCCGAAATG	775
editing	plastid	<i>ndhG</i>	F2	TTTGCTGGACCAATACATG	R2	AGCCACAGAAATTCGACCTAT	517
editing	plastid	<i>petL</i>	F	AAAAAAACATA TTTTATGAGTCCCTTCATG	R	GACCAA TAAACA GA A CTGAGGTTATAG	101
editing	plastid	<i>psbE</i>	F2	ATGCTGGAA GCACAGGAGA	R	CTAAAA CGA TCTACTAAA TTCA TCGAG	252
editing	plastid	<i>psbF</i>	F	GATAGGACCTATCCAAATTTTACAG	R	CGTTGGA TGAACTGCA TTGC	107
editing	plastid	<i>psbZ</i>	F2	TGCTTTCAAATGGCGATTTT	R2	TCAAGAGA TAAGAGAA TTAAGGATACC	181
editing	plastid	<i>rpl23</i>	F2	ATGGATGGAATCAAAATATGC	R2	TTAAGTCTTTTCTTCTAAAGAGG	282
editing	plastid	<i>rpoA</i>	F2	CTCGGACACTACATGGGAAGTGTG	R3	TTCTACGTGAAAATGTTCAA TTTTGATAAG	923
editing	plastid	<i>rpoB</i>	65F	GGTTTTA TCGGTTTATGATCA GGG	792R	CGGCGA CCAA TCCTCTCTAA TTCA C	727
editing	plastid	<i>rpoB</i>	F2	TTCA GGTA TCGACTTCAA AAGAAAC	R2	GGATACTGGGTTCAAATA CCC	472
editing	plastid	<i>rpoC1</i>	F	CCTACTTCTTACGATTA CGA GGTT	862R	ATGGGTCTCA A CTGGGGA GG	374
editing	plastid	<i>rps12-intron</i>	F	ATGCCAACCA TTAACA AACTT	R	TGATGGATTTGACCAA TGGA AAC	295
editing	plastid	<i>rps14</i>	F	TTGATTTA TAGGGA GA A GA GAGGC	R	CCTGGCAA CAAACA TGCCTGAAC	266
splicing	plastid	<i>ndhA</i>	F1	ATGCAACAGCATCCAAACTATAA	R1	GAGTTGAAA GAAGTGTTAATAA	1061
splicing	plastid	<i>petB</i>	F1	GCATTTGATA TTTCCGGAATATGA	R1	GACCCGCCAA TAA CCAA TTT	464
splicing	plastid	<i>petD</i>	F1	CTATGAA GAGATAA TGGA TTATGGG	R1	AGACCTAAAGTTA GAGA TTTATCAA	496
splicing	plastid	<i>rpl2</i>	F1	TGGCGATACATTTATACAAACTT	R1	TACGCGACGAA GAATCAAAGTCT	816
splicing	plastid	<i>rpl16</i>	F1	TCCTTTGATA TAA TTGCTATGCT	R1	AATAATGAA TTGGGTTTATAGG	417
splicing	plastid	<i>ycf3</i>	F1	TCGCGTATAAATGGAAA TTTTAT	R1	GAA GCGCCTCGTATCG	492
editing	mitochondrion	<i>atp1</i>	F2	AGAGCTGCGGA ACTAACGAA TCTA	R5	TGATTTGAGGATTCCTAATGTGATG	1537
editing	mitochondrion	<i>atp4 (orf25)</i>	F1	TGAGATTTGAGTATTA CGAA TATGG	R1	CCCCCGA ACCA TTCTTA	563
editing	mitochondrion	<i>atp6</i>	F2	AA TGCGA CGAA TCTTTTGTGTA	R2	CCAA CGGA ACCA CTAAA TCTT	994
editing	mitochondrion	<i>atp8 (orfB)</i>	F1	GCCTCAACTGGATAAAT	R1	TGTA CAA CATGGATTA GC	453
editing	mitochondrion	<i>atp9</i>	F1	ATGACAAAGCGTGAGTATAAT	R1	TCAGAA TACGAATAAGATCAAA	258
editing	mitochondrion	<i>ccmB (ccb206)</i>	F3	CA GCCTTGAAGTGAATGAAAT	R2	TTAATCTTGTA AACTAATCGAGACC	661
editing	mitochondrion	<i>ccmC (ccb256)</i>	F2	CTACGCGCAAA TTCTCATTTGG	R2	GAGCGAGTGA ACTAAGTTTGGTA	680
editing	mitochondrion	<i>ccmFc (ccb452)</i>	F2	GGTCCAACTACATAACTTTTCTT	R2	ATTA TGA ACTCCA CGGAACTTTCT	1350
editing	mitochondrion	<i>ccmFn-1 (ccb382)</i>	F5	ATTTCCGGGTCTTTCGTGCA TTCACTTA	R5	CGCCGCCGAA GCA GCATGAG	1050
editing	mitochondrion	<i>ccmFn-2 (ccb203)</i>	F2	TGGACACGGGGAGGGAGCAG	R2	CATAACATA ACGGGGGGGGTTGC	601
editing	mitochondrion	<i>cob</i>	F1	GACTATAAGGAACCAACGAT	R1	ATCAGTCTCA TCGGTGTAAG	1171
editing	mitochondrion	<i>cox2</i>	F2	ATGATTTGTTCTAAAATGGTTATT	R1	GTTGGGGGATTAATGATTTGG	770
editing	mitochondrion	<i>matR</i>	F3	GGGCCCTTACTCCGTCCACAC	R2	GTCTTGA CCGGGTCCGAGCTTCC	1651
editing	mitochondrion	<i>mitB</i>	F2	CAC TTTT GCTTGAATTA CTATT	R2	ATTGATAGTTACTTTGCCAGGTTT	852
editing	mitochondrion	<i>nad1</i>	F2	CCA GCTGAAA ATCTTGGAA TAA T	R1	AAA GGTGACTAAA GA CCA GAAAC	947
editing	mitochondrion	<i>nad2</i>	F4	GACGTAACGTAAGTCACTCAGTG	R4	ACGGCCTACCTTCTTTGAA	1578
editing	mitochondrion	<i>nad3</i>	F2	AGCAAGGAGCGA GAAACA AAGT	R2	CCCCATTTTGTGCCATTC	414
editing	mitochondrion	<i>nad4</i>	F4	ATGTTAGAACA TTTCTGTGAA TG	R1	TTTGCCATGTTGCACTAAGTTACT	1479
editing	mitochondrion	<i>nad4L</i>	F1	AA TATTTCA CATTTCTA TGA T	R1	AA TGCTA TTA TAAAT TCTACAG	281
editing	mitochondrion	<i>nad5</i>	F9	ATCAGAAAGGAA GCGCTATAA TGACCA	R8	TAAAAACTACTACTATCAA AATGAAAG	1917
editing	mitochondrion	<i>nad9</i>	F1	TGGATAA CCAA TTCA TTTTCAA	R1	CGTCGCTACGCTGTCC	567
editing	mitochondrion	<i>orf114</i>	F1	CAACCGCGATTTGGATG	R2	AA TGGA AAA GGAACA CCGAGTAG	381
editing	mitochondrion	<i>orf240A</i>	F1	GAGGGGAA GGTGTGTCATAAAT	R1	CGCTGGCGAAAGATA CGAA	636
editing	mitochondrion	<i>rpl2</i>	F2	AGCACTTCTCTA TGGGCA TTGTA	R2	TGCTGCTTCTAATTTGA TGACTGG	811
editing	mitochondrion	<i>rpl5</i>	F1	ATGTTTCCACTCAA TTTTCA TTAC	R1	GAGTTTCCCTCATCTTTT	553
editing	mitochondrion	<i>rps3</i>	F1	AA TCCGATTTCCGGTAAGACTT	R1	CGTTTCCGATATA GCA CGTC	1643
editing	mitochondrion	<i>rps4</i>	F1	TCCCA TTAAGATTTCAA AACTGTGCTG	R1	TTA TATGTTTTGGCCAGTCCGTTTCT	1047
editing	mitochondrion	<i>rps12</i>	F1	TCAA TTGA TTGCTA TGGTAGAGA	R1	TCATA TCGATTTGGGTTTCTTGC	364
editing	mitochondrion	<i>rps14</i>	F1	ATGTCGTA GAA GCA AAA TAGTAGA	R1	TTA TGCCCA TCAA GAA ACC	285
editing	mitochondrion	<i>nad5-1st-intron</i>	F1	GTGGGCGAGGGCTCGTAGTACC	R1	CGGTGGGCTA TCGAACACAGAT	406
editing	mitochondrion	<i>nad7-2nd-intron</i>	F1	CTCCGCCGGTGA CTAAGAAAG	R1	AGCGTGTCTTGGGCCA TCA TAG	343
editing	mitochondrion	<i>cox3</i>	F1	GAGGCA TCTTATCA TTTG	tr-R1	TTCCGGTCA TTTCTTGTGTAAC	994
editing	mitochondrion	<i>nad6</i>	ld-F1	AA GGGCTTGGAA GAA GAA AATG	R2	CTTTTCA CCTTAGTA GTCCATA TGC	653
editing	mitochondrion	<i>nad7</i>	ld-F1	GGAGATGCA TTTCTGGTCAAGTG	R2	ATCCACTCTCAA ACACAATA	1258
editing	mitochondrion	<i>rps7</i>	ld-F1	ACAAACTCGACTAAA GAAGAGGT	R1	ATCTGAAATGCGCGAAACTT	692
editing	mitochondrion	<i>rpl16</i>	F1	AAAA TCGATTATGCTCTG	tr-R1	CCATACATA TCGAGGGCTTATCA	943

Primers combination for multiplex PCR

Primer combination	organelle	gene	size
1	M	nad1	947
		ccb256	680
		orf114	381
2	M	cox3	994
		rpl2	811
		orf240A	636
3	M	rps7	692
		atp4	563
		nad3	414
4	M	nad9	567
		rps12	364
		rps14	285
5	M	cob	1171
		atp6	994
		ccb206	661
6	M	rps4	1047
		nad6	653
7	M	rpl5	553
		nad5-1i	406
		nad7-2ndi	343
		atp9	258
8	M	nad7	1258
		rpl16	943
		cox2	770
9	M	rps3	1643
		nad4	1479
10	M	nad2	1578
		ccb452	1350
11	M+P	nad4L	281
		psbF	107
12	M	atp1	1537
		atp8	453
13	M	ccb382	1050
		ccb203	601
14	P	ndhD	1637
		ndhB	1504
15	P	accD	1638
		matK	1501
16	P	rpl2	816
		rpoB65	727
		ndhG	517
17	P	atpF	520
		rpl16	417
		rps14	266
18	P	clpP	581
		rpoB	472
		psbE	252
19	P	rpl23	282
		psbZ	181
		petL	101
20	P	ndhA	1061
		rpoA	923
		petD	496
21	P	ndhF	775
		ycf3	492
		rpoC1	374
22	P	petB	464
		rps12-i	295
23	M	matR	1651
24	M	nad5	1917
25	M	mttB	852